

## Transport trolley for loading carts TW



Universally usable trolley for safe transportation of MMM loading carts with special guide rolls.

Trolley with running gear, transport plate, handle bar and locking device made of ni-chrome steel. Four swivel castors made of special rubber, with two of them blockable in direction and two totally blockable rolls.

Running gear is adjustable to loading height.

Trolley equipped with safety device for automatic docking of the trolley to the MMM loading and unloading system for sterilisers; easy to unlock manually.

Transport safety device for loading cart, also easy to unlock manually.

Transport trolley	TW 666 EFH 625	TW 666 EFH 850	TW 669 EFH 625	TW 669 EFH 850	TW 6612 EFH 625	TW 6612 EFH 850
Application according to chamber size and loading height in series:	Euro-Selectomat, Selectomat HP, Selectomat Compact/ Securex, Vakulab HP, Vakulab Compact, Seculab, Sterivap		Euro-Selectomat, Selectomat HP, Selectomat Compact/ Securex, Vakulab HP, Vakulab Compact, Ventilab, Seculab, Sterivap		Euro-Selectomat, Selectomat HP, Selectomat Compact/ Securex, Vakulab HP, Sterivap	
Appropriate for	BW 666 Spezial BW 666 Universal		BW 669 Spezial BW 669 Universal		BW 6612 Universal 2 x BW 666	
Outer dimensions with handle bar in mm (Height x width x depth)	875x555x850	1100x555x850	875x555x1070	1100x555x1070	875x555x1290	1110x555x1450
Empty weight in kg approx.	31	32	35	36	37	37
Capacity in sterilisation units STU	4	4	6	6	8	8
Material number	0120573	0109369	0100645	0109518	0109624	0143808



Transport trolley	TW 966 EFH 390	TW 966 EFH 390 Super	TW 969 EFH 390	TW 969 EFH 390 Super	TW 9612 EFH 390	TW 9612 EFH 390 Super
Application according to chamber size and loading height in series:	Vakulab HP		Euro-Selectomat, Selectomat HP, Vakulab HP, Ventilab		Euro-Selectomat, Selectomat HP, Vakulab HP, Seculab, Securex	
Appropriate for	BW 966 Konvect		BW 969 Konvect BW 969 Universal		BW 9612 Konvect BW 9612 Universal 2 x BW 966	
Outer dimensions with bar in (Height x width x depth)	950x555x1010	1010x700x1200	950x555x1170	1010x700x1400	950x555x1440	1010x700x1700
Empty weight in kg approx.	28	42	33	37	53	45
Capacity in sterilisation units STU	6	6	9	9	12	12
Material number	0142957	0143543	0129516	0143905	0143987	0118550

**Lifting trolley** for the transportation of loading carts.

Trolley with smooth transport plate, panelling of telescope column, easy to grab handles and running gear made of ni-chrome steel with ground surfaces.

Automatic lifting mechanism with central telescope column for adjusting the ergonomically best loading and unloading height. Electromechanical lifting spindle drive with rechargeable NiCad batteries for independent power supply.

Foot-switch for automatic, continuous height adjustment.



Docking device for coupling and un-coupling of the trolley to the loading and unloading system of MMM sterilisers. Transport locking of loading carts is easy to unlock.

Running gear consists of four guide rolls of 125 mm diameter with synthetic tire equipment. Two rolls are blockable.

Upper transport plate can be adjusted in height with lower and upper limitation, to be set continuously by limit switch.

Transport trolley	TW-HUB 669 EFH 850/625	TW-HUB 6612 EFH 850/625
Application according to chamber size and loading height in the following series:	Euro-Selectomat, Selectomat HP, Selectomat Compact/ Securex, Vakulab HP	Selectomat Compact/ Securex,
Appropriate for	BW 666/669 Spezial BW 666/669 Universal	BW 6612 Universal 2 x BW 666
Outer dimensions with handle in mm (Height x width x depth)	1100x555x1171	1100x555x1560
Empty weight in kg approx.	112	125
Capacity in sterilisation unit StU	6	8
Material number	0148168	0148557